

**Ultra Morphological and Chemical  
Characterization of an Acid-Base Resistance Zone  
Adjacent to the hybrid layer at Different Time  
Intervals**

*Thesis*

*Submitted to the Faculty  
of Oral and Dental Medicine  
Cairo University*

*In partial fulfillment of Doctor's Degree  
"Operative dentistry"*

**BY**  
**Anwaar Attia Saber**  
**B.D.S. 2003 (AIN SHAMS UNIVERSITY)**  
**M.S. 2009 (CAIRO UNIVERSITY)**

**2012**

# ***Supervisors***

***DR.HEBA SALAH EL-DIN HAMZA***

**Professor of Operative Dentistry  
Faculty of Oral and Dental Medicine  
Cairo University**

***DR . AHMED ZOHEIR EL-HOSHY***

**Assistance Professor of Operative Dentistry  
Faculty of Oral and Dental Medicine  
Cairo University**



*To my parents*

*For their endless support and encouragement*

*To my husband, sister and kids*

*For their endless love and patience*

First, grace and foremost thanks to ***GOD*** for blessing this work until it has reached its end, as a little part of his generous help throughout our life.

I would like to express my sincere gratitude and respect to ***DR. HEBA SALAH EL-DIN HAMZA***, Professor of Operative Dentistry, Faculty of oral and Dental Medicine, Cairo University, for her continuous guidance and supervision, her kind encouragement and support throughout the entire period of the study.

I am deeply grateful to ***DR. AHMED ZOHEIR EL-HOSHY***, Assistance Professor of Operative Dentistry, Cairo University, for his continuous supervision, great help and effort.

I would like also to acknowledge all the staff members of the Operative Dentistry Department, Faculty of Oral and Dental Medicine, Cairo University.

## List of Contents

• List of Tables .....	i
• List of Figures.....	ii
• Introduction.....	1
• Review of literature.....	4
• Aim of the study .....	60
• Materials and Methods .....	61
• Results .....	94
• Discussion.....	140
• Summary and Conclusions.....	148
• References.....	150
• Arabic Summary	

## List of Tables

	Page
Table (1): Variables of the study .....	61
Table (2): Interaction between variables for the grey scale surface roughness and test .....	62
Table (3): Descriptive statistics and test of significance for the effect of group on Ca weight % within each time .....	129
Table (4): Descriptive statistics and test of significance for the effect of time on Ca weight % within each group .....	130
Table (5): Descriptive statistics and test of significance for the effect of group on P weight % within each time .....	131
Table (6): Descriptive statistics and test of significance for the effect of time on P weight % within each group .....	132
Table (7): Descriptive statistics and test of significance for the effect of group on $\mu$ -Tensile Stress (MPa) within each time .....	136
Table (8): Descriptive statistics and test of significance for the effect of time on $\mu$ -Tensile Stress (MPa) within each group .....	137
Table (9): Descriptive statistics and test of significance for the effect of group on flexural bond strength (MPa) within each time .....	138
Table (10): Descriptive statistics and test of significance for the effect of time on flexural bond strength (MPa) within each group .....	139

## List of Figures

	Page
Figure (1): CLEARFIL™ Protect Bond (Kuraray Medical, Tokyo, Japan) .....	67
Figure (2): CLEARFIL™ SE Bond (Kuraray Medical, Tokyo, Japan.....	67
Figure (3): Adper™ Single Bond Adhesive (EB; 3M ESPE, St. Paul, MN, USA) .....	67
Figure (4): Filtek™ Z350 XT Universal restorative(EB; 3M ESPE, St. Paul, MN, USA). ....	68
Figure (5): Bluephase <sup>R</sup> C8 ( Ivoclar Vivadent AG, Schaan/liechtensein, Austria .....	68
Figure(6): The buffered demineralizing solution .....	68
Figure (7): Human molars were collected and thoroughly examined to be free from caries, hypocalcification or any restoration. ....	77
Figure (8): The coronal portion of each tooth was removed using a low speed diamond saw under water cooling (Isomet 1000, Buehler Ltd., Lake Bluff, IL, USA) to expose a flat midcoronal dentin surface.....	77
Figure(9): The coronal portion of each tooth was removed using .....	77
Figure(10): Silicon Carbide abrasive paper .....	78



Figure (11): Putting the primer of clearfil protect bond.....	78
Figure (12): Putting the adhesive of clearfil protect bond .....	78
Figure (13): Applying Clearfil Protect Bond primer .....	79
Figure (14): Gently air blow the primer of Clearfil protect Bond.....	79
Figure(15): Applying Clearfil Protect Bond adhesive .....	79
Figure (16): Gently air blow the adhesive of Clearfil protect Bond.....	80
Figure (17): The specimen was light cured for 10 seconds using the light curing system bluephase <sup>R</sup> C8 .....	80
Figure(18): The primer and adhesive of Clearfil SE bond.....	80
Figure(19): Putting the primer of Clearfil SE bond.....	81
Figure(20): Putting the adhesive of Clearfil SE Bond .....	81
Figure(21): Applying Clearfil SE Bond primer .....	81
Figure (22): Gently air blow the primer of Clearfil SE Bond.....	82
Figure(23): Applying Clearfil SE Bond.....	82
Figure(24): Gently air blow the adhesive of Clearfil SE Bond.....	82
Figure(25): The specimen was light cured for 10 seconds using the light curing system bluephase <sup>R</sup> C8 .....	83

Figure(26): Applying scotchbond etchant 35% Phosphoric acid .....	83
Figure(27): Washing the etchant with water .....	83
Figure(28): Gently air blow the excess water .....	84
Figure (29): Applying Single Bond two coats.....	84
Figure (30): Gently air blow Single Bond adhesive .....	84
Figure (31): The specimen was light cured for 10 seconds using the light curing system bluephase <sup>R</sup> C8 .....	85
Figure(32): Building up the composite .....	85
Figure(33): Using auto matrix for building up the composite .....	85
Figure(34): Filtek Z-350 was built up in one increment approximately 1-mm thick for SEM testing.....	86
Figure(35): 1mm <sup>2</sup> increment was light cured for 40 seconds using the light curing system bluephase <sup>R</sup> C8 .....	86
Figure(36): The fifth increment was applied .....	86
Figure(37): Putting the matrix down and curing.....	87
Figure(38): Light curing the fifth increment .....	87
Figure(39): Clearfil Protect Bond specimen for Microtensile strength test.....	87
Figure(40): Clearfil SE Bond specimen for Microtensile strength test.....	88

Figure (41) Single Bond specimen for Microtensile strength test:.....	88
Figure(42): The cutting machine with diamond disc mounted .....	88
Figure(43): The tooth mounted on the cutting machine for separating into two halves.....	89
Figure(44): Diamond disc mounted with separating wheel spacer .....	89
Figure(45): The prepared specimen stored in buffering solution for 90 minutes .....	89
Figure(46): Gold sputtering for the SEM specimens .....	90
Figure(47): The Scanning Electron Microscope for 2 samples. Using SEM Model Quanta FEG attached with EDAX Unit .....	90
Figure (48): The tooth after cutting into beam sticks and before separating the beams. Some central beams which had the pulp chamber underneath were freely dropped before separating the beams.....	90
Figure (49): A beam after separation consisting of a composite and dentin parts .....	91
Figure(50): Beam width and breadth measurement before $\mu$ TBS testing .....	91
Figure (51): The attachment jig (front view) .....	91
Figure (52): Rear surface of the attachment jig: the arrow points to the firm plastic sheet.....	92
Figure (53): Diagram showing the components of the attachment jig:A: is the fixed part; while B: is the moving component.....	92

Figure (54): The attachment jig face showing a glued beam in the orientation groove during testing.....	92
Figure(55): Lloyd universal testing machine .....	93
Figure (56): Flexural Bond strength (Three point bend test set up). .....	93
Figure (57): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (1000x).....	98
Figure(58): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (2000x).....	98
Figure(59): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (3500x) .....	98
Figure (60): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (1000x).....	99
Figure (61): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (2000x).....	99
Figure (62): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (3500x) .....	99
Figure (63): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (1000x).....	100
Figure(64): SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (2000x).....	100

<b>Figure (65):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after 24 hours subjected to acid-base challenge (3500x).....	100
<b>Figure (66):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after three months subjected to acid-base challenge (500x).....	101
<b>Figure(67):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after three months subjected to acid-base challenge (1000x) .....	101
<b>Figure(68):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after three months subjected to acid-base challenge (7000x) .....	101
<b>Figure(69):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after three months subjected to acid-base challenge (1000x) .....	102
<b>Figure(70):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after three months subjected to acid-base challenge (3500x) .....	102
<b>Figure(71):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after three months subjected to acid-base challenge (7000x) .....	102
<b>Figure(72):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after six months subjected to acid-base challenge (1000x) .....	103
<b>Figure(73):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after six months subjected to acid-base challenge (3500x) .....	103

<b>Figure(74):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after six months subjected to acid-base challenge (7000x) .....	<b>103</b>
<b>Figure(75):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after six months subjected to acid-base challenge (1000x), .....	<b>104</b>
<b>Figure(76):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after six months subjected to acid-base challenge (3500x) .....	<b>104</b>
<b>Figure(77):</b> SEM observations of the adhesive-dentin interfaces for Clearfil Protect Bond after six months subjected to acid-base challenge (5000x), .....	<b>104</b>
<b>Figure(78):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after 24 hours subjected to acid-base challenge (1000x) .....	<b>105</b>
<b>Figure(79):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after 24 hours subjected to acid-base challenge (2000x) .....	<b>105</b>
<b>Figure(80):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after 24 hours subjected to acid-base challenge (3500x).....	<b>105</b>
<b>Figure(81):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after 24 hours subjected to acid-base challenge (1000x).....	<b>106</b>
<b>Figure(82):</b> SEM observations of the adhesive-dentin interfaces for	

Clearfil SE Bond after 24 hours subjected to acid-base challenge (2000x) .....	106
<b>Figure(83):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after 24 hours subjected to acid-base challenge (3500x) .....	106
<b>Figure(84):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after three months subjected to acid-base challenge (1000x).....	107
<b>Figure(85):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond Bond after three months subjected to acid-base challenge (2000x) .....	107
<b>Figure(86):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after three months subjected to acid-base challenge (3500x) .....	107
<b>Figure(87):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after three months subjected to acid-base challenge (1000x).....	108
<b>Figure(88):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after three months subjected to acid-base challenge (2000x) .....	108
<b>Figure(89):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after three months subjected to acid-base challenge (3500x) .....	108
<b>Figure(90):</b> SEM observations of the adhesive-dentin interfaces for Clearfil SE Bond after six months subjected to acid-base challenge (1000x).....	109